




















LISTADO DE CIRCUITOS DE PACIENTE

No. DE PARTE	DESCRIPCIÓN	IMAGEN
CIRCUITOS DE ANESTESIA TIPO BAIN		
8570073	ADU BAIN CIRCUIT, COAXIAL MAPLESON D CIRCUIT (1.5 M)	
8570074	AESTIVA® BAIN CIRCUIT, COAXIAL MAPLESON D CIRCUIT (1.5 M), WITH GAS SAMPLING ELBOW	
8570076	MODIFIED JACKSON REES CIRCUIT KIT, MODIFIED MAPLESON F CIRCUIT WITH TAIL VALVE, EXPIRATORY LIMB 0.25 M LONG, FRESH GAS HOSE 1.8 M LONG. INCLUDES: <ul style="list-style-type: none"> • 0.5 LITER LATEX-FREE MANUAL BAG NOTE: FITS ALL DATEX-OHMEDA ANESTHESIA MACHINES	
CIRCUITOS DE ANESTESIA TIPO UNIT		
M1007452	UNI-CIRCUIT KIT, PEDIATRIC, 1.5 M/60 IN., INNER TUBE DIAMETER 10 MM, OUTER TUBE DIAMETER 22 MM, CUFF DIAMETER: 22 MM, DEAD SPACE: 17 ML, TOTAL VOLUME: 680 ML, RESISTANCE TO FLOW: EXPIRED 4.2 CM H2O AT 60 L/MIN. FLOW, INSPIRED 13.9 CM H2O AT 60 L/MIN. FLOW, EXPIRED 1.14 CM H2O AT 30 L/MIN. FLOW, INSPIRED 3.51 CM H2O AT 30 L/MIN. FLOW. INCLUDES: <ul style="list-style-type: none"> • GAS SAMPLING ELBOW • 1 LITER LATEX-FREE MANUAL BAG 	
8570075	UNI-CIRCUIT, ADULT, 1.5 M/60 IN. INNER TUBE DIAMETER 15 MM, OUTER TUBE DIAMETER 25 MM, DEAD SPACE: 17 ML, TOTAL VOLUME: 1100 ML, RESISTANCE TO FLOW: EXPIRED 1.9 CM H2O, INSPIRED 2.8 CM H2O AT 60 L/MIN. FLOW. INCLUDES: <ul style="list-style-type: none"> • GAS SAMPLING ELBOW 	
CIRCUITOS DE ANESTESIA PARA MÁQUINA ADU		
M1035988	ADU MANUAL VENTILATION KIT. INCLUDES: <ul style="list-style-type: none"> • MANUAL VENTILATION TUBE (1.5 M/60 IN.) WITH QUICK CONNECTOR • VENTILATION TUBE (1.0 M/40 IN.) NOTE: REPLACES REF 8570183	
8570085	ADU MANUAL VENTILATION KIT. INCLUDES: <ul style="list-style-type: none"> • MANUAL BAG TUBE (1.5 M/60 IN.) WITH QUICK CONNECTOR • 2 LITER LATEX-FREE MANUAL BAG 	





LISTADO DE CIRCUITOS DE PACIENTE

No. DE PARTE	DESCRIPCIÓN	IMAGEN
M1019528	<p>ADU MANUAL VENTILATION KIT, ADULT, 1.5 M/60 IN. INCLUDES:</p> <ul style="list-style-type: none"> • MANUAL VENTILATION TUBE WITH QUICK CONNECTOR • 2 LITER LATEX-FREE MANUAL BAG <p>NOTE: REPLACES REF 8570085</p>	
M1019508	<p>ADU PATIENT CIRCUIT KIT, PEDIATRIC, 1.5 M/60 IN., TUBING DIAMETER 15 MM WITH 22 MM CUFFS. INCLUDES:</p> <ul style="list-style-type: none"> • STANDARD ELBOW • MANUAL VENTILATION TUBE (1.5 M/60 IN.) WITH QUICK CONNECTOR • VENTILATION TUBE (1.0 M/40 IN.) • 1 LITER LATEX-FREE MANUAL BAG <p>NOTE: REPLACES REF M1010865</p>	
M1035985	<p>ADU PATIENT CIRCUIT KIT, PEDIATRIC, 1.5 M/60 IN. INCLUDES:</p> <ul style="list-style-type: none"> • STANDARD ELBOW • MANUAL VENTILATION TUBE (1.5 M/60 IN.) WITH QUICK CONNECTOR • VENTILATION TUBE (1.0 M/40 IN.) • 2 LITER LATEX-FREE MANUAL BAG <p>NOTE: REPLACES REF 8570166</p>	
M1019529	<p>ADU PATIENT CIRCUIT KIT, ADULT, 1.5 M/60 IN. INCLUDES:</p> <ul style="list-style-type: none"> • STANDARD ELBOW • MANUAL VENTILATION TUBE (1.5 M/60 IN.) WITH QUICK CONNECTOR • VENTILATION TUBE (1.0 M/40 IN.) • 2 LITER LATEX-FREE MANUAL BAG <p>NOTE: REPLACES REF 8570086</p>	
887420	<p>HYTREL MANUAL VENTILATION TUBE, 1.4 M/56 IN., CONNECTORS: 22 MM M, QUICK CONNECTOR TO THE BELLOWS BLOCK, FOR AS/3, S/5 ADU</p>	





LISTADO DE CIRCUITOS DE PACIENTE

No. DE PARTE	DESCRIPCIÓN	IMAGEN
CIRCUITOS DE ANESTESIA		
8002930	PATIENT CIRCUIT KIT, ADULT, 1.0 M/40 IN., TUBING DIAMETER 22 MM INCLUDES: • GAS SAMPLING ELBOW • 3 LITER LATEX-FREE MANUAL BAG • ANESTHESIA GAS SAMPLING LINE (3.0 M/10 FT.).	
M1019547	PATIENT CIRCUIT KIT, ADULT, 1.5 M/60 IN. INCLUDES: • GAS SAMPLING ELBOW • 2 LITER LATEX-FREE MANUAL BAG NOTE: REPLACES REF 8570118	
8570114	PATIENT CIRCUIT, PEDIATRIC, 1.5 M/60 IN., TUBING DIAMETER 15 MM INCLUDES: • HMEF 500	
M1019539	PATIENT CIRCUIT, ADULT 1.8 M/72 IN. INCLUDES: • STANDARD ELBOW • HMEF 1000/S NOTE: REPLACES REF 8570103	
8570096	PATIENT CIRCUIT KIT, ADULT, 1.8 M/72 IN., TUBING DIAMETER 22 MM INCLUDES: • GAS SAMPLING ELBOW • 2 VENTILATION TUBES (1.0 M/40 IN.) • 2 LITER LATEX-FREE MANUAL BAG	
M1019534	PATIENT CIRCUIT KIT, ADULT, 1.8 M/72 IN. INCLUDES: • GAS SAMPLING ELBOW • MANUAL VENTILATION TUBE (1.0 M/40 IN.) • VENTILATION TUBE (1.0 M/40 IN.) • 2 LITER LATEX-FREE MANUAL BAG NOTE: REPLACES REF 8570096	
8004457	PATIENT CIRCUIT KIT, ADULT, FOR ABS WITH BAG ARM, 1.8 M/72 IN., TUBING DIAMETER 22 MM. INCLUDES: • GAS SAMPLING ELBOW • MANUAL VENTILATION TUBE (0.9 M/36 IN.) WITH ABS BAG ARM CONNECTOR • 3 LITER LATEX-FREE MANUAL BAG • ANESTHESIA GAS SAMPLING LINE (3.0 M/10 FT.)	







LISTADO DE CIRCUITOS DE PACIENTE

No. DE PARTE	DESCRIPCIÓN	IMAGEN
8004458	PATIENT CIRCUIT KIT, ADULT, 1.8 M/72 IN., TUBING DIAMETER 22 MM. INCLUDES: <ul style="list-style-type: none"> • GAS SAMPLING ELBOW • MANUAL VENTILATION TUBE (1.8 M/72IN.) • 2 LITER LATEX-FREE MANUAL BAG • STRAIGHT CONNECTOR BETWEEN BAG AND TUBE • ANESTHESIA GAS SAMPLING LINE (3.0 M/10 FT.) 	
8004461	SCAVENGING TUBE, SPECIFIC FOR ADVANCED BREATHING SYSTEM (ABS), TUBING DIAMETER 30 MM, CONNECTORS: 30 MM F, LENGTH 5.0 M/197 IN.	
M1093632	ABS KIT ROW WITH BAG SUPPORT ARM DISPOSABLE, INCLUDES: <ul style="list-style-type: none"> • 72 IN (1.8 M) CORRUGATED HOSES • 2 L LATEX-FREE BREATHING BAG • PARALLEL WYE • ELBOW WITH CAPPED GAS SAMPLING PORT • 36 IN (0.9 M) MANUAL VENTILATION HOSE • ABS BAG ARM CONNECTOR • 7 IN (18 CM) EXHAUST LINE WITH COLDER FITTING FOR GAS RETURN OR SCAVENGING • 40 IN (1.0 M) EXHAUST LINE WITH COLDER FITTING FOR GAS RETURN OR SCAVENGING • 10 FT (3.0 M) CO2 LINE 	
M1093633	ABS KIT ROW FOR BAG ON A TUBE, DISPOSABLE, INCLUDES: <ul style="list-style-type: none"> • 72 IN (1.8 M) CORRUGATED HOSES • 2 L LATEX-FREE BREATHING BAG • PARALLEL WYE • ELBOW WITH CAPPED GAS SAMPLING PORT • 72 IN (1.8 M) MANUAL VENTILATION HOSE • 7 IN (18 CM) EXHAUST LINE WITH COLDER FITTING FOR GAS RETURN OR SCAVENGING • 40 IN (1.0 M) EXHAUST LINE WITH COLDER FITTING FOR GAS RETURN OR SCAVENGING • 10 FT (3.0 M) CO2 LINE 	

LISTADO DE CIRCUITOS DE PACIENTE

No. DE PARTE	DESCRIPCIÓN	IMAGEN
70100	CIRCUITO DE PACIENTE	
CIRCUITOS DE VENTILACIÓN		
M1012148	<p>PATIENT CIRCUIT KIT, PEDIATRIC, 1.5 M/60 IN., TWO LIMB CIRCUIT, FOR ACTIVE HUMIDIFICATION, ICU APPLICATION INCLUDES:</p> <ul style="list-style-type: none"> • DOUBLE SWIVEL ELBOW • Y-PIECE WITH GAS SAMPLING PORT (Y-PIECE IS DETACHABLE FROM CIRCUIT TO INSERT NEBULIZER) • WATER TRAPS IN THE MIDDLE OF BOTH LIMBS • 14" TUBE WITH ELBOW (FOR HUMIDIFIER) 22MM F-22MM F 	 <p>The diagram shows a patient circuit kit for active humidification. It features a central Y-piece with a gas sampling port, two limbs with water traps, and a 14-inch tube with an elbow for the humidifier. The kit includes a double swivel elbow and a 360-degree rotating connector.</p>
M1012152	<p>PATIENT CIRCUIT, PEDIATRIC, 1.5 M/60 IN., TWO LIMB CIRCUIT, FOR PASSIVE HUMIDIFICATION, ICU APPLICATION INCLUDES:</p> <ul style="list-style-type: none"> • DOUBLE SWIVEL ELBOW • Y-PIECE WITH GAS SAMPLING PORT (Y-PIECE IS DETACHABLE FROM CIRCUIT TO INSERT NEBULIZER) 	 <p>The diagram shows a pediatric patient circuit for passive humidification. It includes a Y-piece with a gas sampling port, two limbs, and a double swivel elbow. The circuit is designed for ICU application and features a 360-degree rotating connector.</p>
M1012145	<p>PATIENT CIRCUIT, ADULT, 1.5 M/60 IN. TWO LIMB CIRCUIT, FOR PASSIVE HUMIDIFICATION, ICU APPLICATION INCLUDES:</p> <ul style="list-style-type: none"> • DOUBLE SWIVEL ELBOW • Y-PIECE WITH GAS SAMPLING PORT (Y-PIECE IS DETACHABLE FROM CIRCUIT TO INSERT NEBULIZER) 	 <p>The diagram shows an adult patient circuit for passive humidification. It consists of a Y-piece with a gas sampling port, two limbs, and a double swivel elbow. The circuit is suitable for ICU application and includes a 360-degree rotating connector.</p>
M1012142	<p>PATIENT CIRCUIT KIT, ADULT, 1.5 M/60 IN. TWO LIMB CIRCUIT, FOR ACTIVE HUMIDIFICATION, ICU APPLICATION. INCLUDES:</p> <ul style="list-style-type: none"> • DOUBLE SWIVEL ELBOW • Y-PIECE WITH GAS SAMPLING PORT (Y-PIECE IS DETACHABLE FROM CIRCUIT TO INSERT NEBULIZER) • 14 IN. TUBE WITH ELBOW (FOR HUMIDIFIER) 22 MM F – 22 MM F • WATER TRAPS IN THE MIDDLE OF BOTH LIMBS 	 <p>The diagram shows an adult patient circuit kit for active humidification. It features a central Y-piece with a gas sampling port, two limbs with water traps, and a 14-inch tube with an elbow for the humidifier. The kit includes a double swivel elbow and a 360-degree rotating connector.</p>

LISTADO DE CIRCUITOS DE PACIENTE

No. DE PARTE	DESCRIPCIÓN	IMAGEN
8570101	PATIENT CIRCUIT KIT, ADULT, 2.3 M/92 IN., TUBING DIAMETER 22 MM INCLUDES: <ul style="list-style-type: none"> • VENTILATION TUBE (1.0 M/40 IN.) • Y-PIECE WITH TEMP PORT • 2 WATER TRAPS 	
CIRCUITO DE PACIENTE		
M1019544	PATIENT CIRCUIT, PEDIATRIC, 1.5 M/60 IN. NOTE: REPLACES REF 8570111	
M1019543	PATIENT CIRCUIT, EXPANDABLE, ADULT, 1.5 M/60 IN. NOTE: REPLACES REF 8570109	
8570112	PATIENT CIRCUIT, PEDIATRIC, 1.5 M/60 IN., TUBING DIAMETER 15 MM INCLUDES: <ul style="list-style-type: none"> • 1 WATER TRAP 	
8570111	PATIENT CIRCUIT, PEDIATRIC, 1.5 M/60 IN., TUBING DIAMETER 15 MM, DEAD SPACE: 1.5 ML, TOTAL VOLUME: 680 ML	
CIRCUITOS DE ANESTESIA PARA USO EN MRI		
M1040965	PATIENT CIRCUIT, ADULT, 6.0 M/240 IN., TWO LIMB CIRCUIT FOR USE IN MRI ENVIRONMENTS. INCLUDES: <ul style="list-style-type: none"> • GAS SAMPLING ELBOW 	
CIRCUITOS DE VENTILACIÓN PARA SISTEMA DE ENTREGA APTAÉR HELIOX		
M1004517	HELIOX PATIENT CIRCUIT, ADULT, 1.5 M/60 IN. INCLUDES: <ul style="list-style-type: none"> • EXHALATION VALVE PERMANENTLY ATTACHED • PRESSURE SENSING LINE (1.8 M/72 IN.) 	